

CLINICAL LAB DEVICES

HC-B018

Microplate Washei

Application

- Hospitals and relevant departments of the enzyme standard testing enzymes panels, training board flushing work
- Chemical, biochemical, etc in the aspects of trace analysis corresponding cleaning

Features

- Liquid crystal display
- Convenient setting and operation mode
- Save preset programs in the memory
- Low noise working mode

Main Technical Parameters

Item	Contents
Washing Head	8 or 12 channels
Applicable Microplate	96 or 48 wells microplate
Washing Times Option	1-20,presettable
Washing Rows Option	1-12,presettable
Soak Time	1~3600s
Residual Volume	<2µL/well



Guangzhou Happycare Electronics Co., Ltd.

Phone: +86-15820241249

Email : sale001@happycaregroup.com Website : www.happycaregroup.com

Address: Room 2202–2203, Dongcheng Square, No. 5, Huangpu East Road, Huangpu District, Guangzhou, Guangdong, China